

Troubleshooting Guide

Problem	Possible Reason	Solution
Poor bacterial growth	Inoculated bacterial sample from a plate or a culture stock stored over a long time period	Always inoculate bacterial cells from a freshly streaked plate and grow with required antibiotic(s).
	Inadequate shaking during incubation	Grow cells with vigorous shaking (e.g. 250 rpm). Adjust a suitable shaking speed according to the angular magnitude of an orbital shaker platform.
Poor cell lysis Low yield of plasmid DNA	Use of excessive amount of bacterial cells harvested from a large or over-grown culture	Up to 5 ml culture for high-copy plasmid. Up to 10 ml culture for low-copy plasmid. When the culture volume is larger than 5 ml, increase the amount of MX1, MX2 and MX3 Buffer.
	Insufficient amount of bacterial cells	Ensure that bacteria have been grown well ($OD_{600} > 1$) after overnight incubation at vigorous shaking.
Low yield of plasmid DNA	Overgrowth of bacteria	Incubate bacterial culture with LB medium and do not incubate for more than 16 hours.
	Plasmid does not propagate	Always inoculate bacterial cells from a freshly streaked plate and grow with required antibiotic(s)
	Inefficient or incomplete DNA elution	Make sure that elution solution is at pH 7-8.5, in full contact with and completely absorbed by the membrane.

Troubleshooting Guide

Problem	Possible Reason	Solution
Low yield of plasmid DNA	Poor cell lysis	Refer to Solution section of problem – “Poor cell lysis”.
	Plasmid is larger than 10-kbp	Use elution solution preheated to 70°C.
Plasmid appears smeared or degraded	Host strain is <i>endA</i> ⁺	Use <i>endA</i> ⁻ strain if possible. Wash with additional 250 µl of MX3 Buffer before proceeding to the WN washing step.
	Overgrowth of bacteria	Incubate bacterial culture with LB medium, and do not incubate for more than 5 minutes.
Genomic DNA contamination in eluate	Lysate improperly prepared	After addition of MX2 Buffer, mix gently to prevent genomic DNA from shearing. Do not incubate for more than 5 minutes.
RNA contamination	Insufficient RNase A activity in MX1 Buffer	Ensure that entire RNase A is added into MX1 Buffer and stored at 4°C. After long-term storage (> 6 months), add RNase A into MX1 Buffer to a conc. of 70µg/ml and store at 4°C.
Plasmid of poor quality	Ethanol in WS buffer is not completely removed	Following WS washing step, discard the flow-through and centrifuge the column at full speed >/= 3 minutes.
	Use of excessive amount of bacterial cells harvested from a large or over-grown culture	Reduce the amount of sample used. Incubate bacterial culture with LB medium, and do not incubate for more than 16 hours.