

Data Sheet (11.08.2010)


mi-20 bp DNA Marker Go

20 bp to 200 bp log scale

ready-to-use

Cat.-No.	Amount	Conc.
mi-M8202	100 rx	0.1 µg/µl

For in vitro use only!



Description

mi-20 bp DNA Marker Go consists of 10 DNA fragments ranging from 20 bp to 200 bp and is supplied in ready-to-use format containing an orange tracking dye. It allows the estimation of DNA fragments on agarose gels generated by PCR or restriction digest.

DNA log scale ladders are well suited for combination if small and large fragments should be analyzed on one gel.

Recommended Load

5 µl DNA ladder per lane

DNA Fragments

20, 40, 60, 80, **100**, 120, 140, 160, 180, 200 bp

Form

DNA fragment mix in 10 mM Tris-HCl (pH 7.5), 10 mM EDTA, 10% (w/v) Glycerol, Orange G

Storage

Store at -20°C.

Short term storage (up to 4 weeks) at 4°C possible.